






LP ITALIANA SPA
























"S" shaped ABS rack, resistant up to 90 °C
 Portaprovette a "S" in ABS, resistente fino a 90 °C
 Portoir à "S" en ABS, résistant jusqu'à 90 °C

	code	dimensions mm	places posti places	tube mm provetta mm tube mm	colour colore couleur
1	183082	225x114x60	32	Ø 21	
	183083	225x114x60	32	Ø 21	
	183084	225x114x60	32	Ø 21	
	183085	225x114x60	32	Ø 21	
	183086	225x114x60	32	Ø 21	
2	184082	225x114x60	50	Ø 17	
	184083	225x114x60	50	Ø 17	
	184084	225x114x60	50	Ø 17	
	184085	225x114x60	50	Ø 17	
	184086	225x114x60	50	Ø 17	
3	185082	180x94x50	50	Ø 13	
	185083	180x94x50	50	Ø 13	
	185084	180x94x50	50	Ø 13	
	185085	180x94x50	50	Ø 13	
	185086	180x94x50	50	Ø 13	



Square hole acetal resin rack, resistant up to 121 °C
 Portaprovette a fori quadri in resina acetalica, resistente fino a 121 °C
 Portoir à trous carrés en résine acétalique, résistant jusqu'à 121 °C

	code	dimensions mm	places posti places	tube mm provetta mm tube mm	colour colore couleur
1	187081	246x101x70	40	Ø 21	
	187082	246x101x70	40	Ø 21	
	187083	246x101x70	40	Ø 21	
	187084	246x101x70	40	Ø 21	
	187085	246x101x70	40	Ø 21	
	187086	246x101x70	40	Ø 21	
2	188081	246x101x70	60	Ø 17	
	188082	246x101x70	60	Ø 17	
	188083	246x101x70	60	Ø 17	
	188084	246x101x70	60	Ø 17	
	188085	246x101x70	60	Ø 17	
	188086	246x101x70	60	Ø 17	
3	189081	246x101x70	90	Ø 13	
	189082	246x101x70	90	Ø 13	
	189083	246x101x70	90	Ø 13	
	189084	246x101x70	90	Ø 13	
	189085	246x101x70	90	Ø 13	
	189086	246x101x70	90	Ø 13	



Square hole polystyrene price saving test tube rack
Portaprovette economico a fori quadrati in polistirolo
 Portoir économique à trous carrés en polystyrène


	code	dimensions mm	places posti places	tube mm provetta mm tube mm
1	187011	246x101x70	40	Ø 21
	188011	246x101x70	60	Ø 17
	189011	246x101x70	90	Ø 13
2	187012	246x101x70	40	Ø 21
	188012	246x101x70	60	Ø 17
	189012	246x101x70	90	Ø 13

White polypropylene rack
Portaprovette in polipropilene bianco
 Portoir en polypropylène blanc

code	dimensions mm	places posti places	tube mm provetta mm tube mm
180042	243x126x64	18	Ø 30



Styro-foam test tube tray
Vassoio in polistirolo espanso per provette
 Portoir pour tube à essai en polystyrène expansé




code	dimensions mm	places posti places	tube mm provetta mm tube mm
290370	207x207x56	100	Ø 12
290371	207x207x56	100	Ø 16

Stackable and non floating ABS rack
with alphanumeric grid

*Portaprovette impilabili non galleggianti in
ABS con griglia alfanumerica*

Portoir empilable et non flottant, en ABS,
grille alphanumérique




	code	dimensions mm	places posti places	tube mm provetta mm tube mm
1	181082	260x130x60	50	Ø 17
2	182082	210x105x50	50	Ø 13

Square hole ABS rack for cuvettes

Portacuvette in ABS a fori quadri

Portoir en ABS à trous carrés pour cuves



code	dimensions mm	places posti places
186082	250x54x37	24