



FOC 200E Connect

Incubador refrigerador

Excelente estabilidad de temperatura y homogeneidad
Puerta transparente interna



Incubador habilitado en nube

- Auto-Tuning en su sistema de termoregulación homogeneidad y estabilidad ($\pm 0,5^{\circ}\text{C}$).
- Selección de temperatura ajustable en el rango de 3 – 50 °C.
- Sistema de enfriado clase A+.
- Equipado con conexión Wi-Fi a la plataforma en nube Velp Ermes.
- Conexión opcional vía cable a PC y software TEMPSoft™.
- MODBUS protocolo comunicación.
- Conexiones internas para alimentación de dispositivos.
- Interface digital intuitiva.



CÓDIGO - INSTRUMENTO

FOC 200E Connect 230V/50-60Hz F10300542



ESPECIFICACIONES TÉCNICAS

PANTALLA DIGITAL	3-dígitos, 0.1 °C resolución
VOLUMEN INTERNO ÚTIL	169 litros
SISTEMA DE TERMOREGULACIÓN ELECTRONICA	Auto-Tuning
RANGO TEMPERATURA	DESDE 3.0 A 50.0°C
ESTABILIDAD DE TEMPERATURA INTERNA	$\pm 0.5^{\circ}\text{C}$
HOMOGENEIDAD DE TEMPERATURA INTERNA	$\pm 0.5^{\circ}\text{C}$
NÚMERO DE BANDEJAS INCLUIDAS	3
NÚMERO DE CONECTORES INTERNOS	2
ERMES CONEXION	via Wi-Fi
PC CONEXION	OpCIONAL via RS485 cable y TEMPSoft™ Software
INTERFACE COMUNICACIÓN	Wi-Fi, 2 x RS485, MODBUS
CONSUMO	400 W
ALIMENTACION	230 V / 50-60 Hz
PESO	46.3 Kg
DIMENSIONES (WXHxD)	540x1300x550 mm

Rev 2 12.2020

05.2017

Velp reserves the right to make technical alteration
Velp does not assume liability for errors in printing, typing or transmission



DESIGNED AND MANUFACTURED IN ITALY



Palex

Constant Improvement

<https://bit.ly/3uPuMI4>